

Region of interest (ROI) selection using vision transformer for automatic analysis using whole slide images

Md Shakhawat Hossain, Galib Muhammad Shahriar, M. M. Mahbubul Syeed, Mohammad Faisal Uddin, Mahady Hasan, Shingla Shivam & SureshAdvani

Department of Computer Science and Engineering, Independent University Bangladesh, Dhaka 1229, Bangladesh

RIoT Research Center, Independent University Bangladesh, Dhaka 1229, Bangladesh

Department of Pathology, SL Raheja Hospital, Mumbai 400016, India

*Email: shakhawat@iub.edu.bd

Extended Abstract:

Goal: The Study aims that the choice of regions of interest (ROI) is a routine step in medical image analysis for all images modalities. also enhancing diagnostic accuracy, compromising inter-observer variability.

Research Area: This research mainly focused by Deep leaning, medical AI and Computational digital image analysis. Here, A whole slide image (WSI) is a digital image of a pathological specimen generated by a WSI scanner and used.

Methodology: There are few methods on finding ROI from WSI mentioned in literature. The methods can broadly be classified in two categories: those that are based on image-features and those that are based on psycho-physical features of pathologists^{13,14} also WSI quality evaluation. Methods of image features make use of statistics or structures.

Results and Analysis: The ViT-based method outperformed CNN models, showing higher accuracy (99% at 20 \times , 97% at 10 \times) and faster processing (less than 15 seconds per WSI). This study investigated the role of image magnification in ROI detection.

References:

[1] Hossain, M. S. et al. "Region of Interest Selection Using Vision Transformer for Automatic Analysis Using Whole Slide Images." *Scientific Reports* 13, 11314 (2023).




Figure 2. Human epidermal growth factor receptor 2 (HER2) status changes depending on the selection of region regardless of the test methods.




Figure 4. Overlap of ROIs for proposed method and pathologists.




Figure 6. ROC for 5-fold cross validation for 10x images.




Figure 7. ROC for 5-fold cross validation for 20x images.

Keywords:

ROI detection, Vision Transformer, Whole Slide Image, Digital Pathology.

